



# Eyewash

First aid in the event of eye accidents



In case of eye accidents caused by caustic liquids or small foreign bodies, rapid reaction is required in order to minimise the risk of permanent damage, and here Plum Eyewash is both efficient and easy to use.

The sterile isotonic saline solution cleanses the eye and washes away foreign bodies or sprays of harmful liquids, thus minimising the amount of damage done. The saline solution's pH approximates to the pH of the eye and caters for a pleasant and effective eye wash with no side effects.

Eyewash is easy to use and ideal for utilisation in industry, production and chemical environments, but should be readily available at any location with an increased risk of eye accidents and damage.

## Eyewash

Prod. no.:	PL 4652
Contents:	500 ml 0.9 % sterile sodium chloride solution
No. of bottles in station:	2 pcs.
Rinse time:	Approx. 10 min.
Extras:	Mirror in Wash station lid

Lotek - Eyewash - 120112



- et selskab i GD-Gruppen

Lotek A/S, Bohrsvej 7, DK-8600 Silkeborg, Tel. +45 70 13 52 00, Fax +45 86 80 32 39, www.lotek.dk, mabe@lotek.dk

